




MAGiC Image Gallery




T1 FLAIR




T2 FLAIR




R1 Map




T1




PD




PSIR



T2



T1 Map



PSIR



Imagination at work

Product may not be available in all countries and regions.
Contact a GE Healthcare Representative for more information.

GE Healthcare
3000 N. Grandview Blvd.
Waukesha, WI 53188
U.S.A.

Date subject to change.

©2015 General Electric Company.

SIGNA, GE, the GE Monogram, and imagination at work
are trademarks of General Electric Company.

Reproduction in any form is forbidden without prior written
permission from GE. Nothing in this material should be used to
diagnose or treat any disease or condition. Readers must
consult a healthcare professional.

MR-0496-10.15-EN-US
JB35506SE


GE Healthcare

Believe

in the power of One & Done imaging.

MAGiC

Imagine what MR can be.



CE marked. Not available for sale in all countries.



Reveal

a new wonderland of MR possibilities.

The secret of MAGiC (MAGnetic resonance image Compilation) lies in its unique ability to make possible multiple image contrasts in a single neuro scan. We call this One & Done imaging, and its benefits are wide-ranging.


It reduces scanning time by as much as 50% allowing you to significantly increase your case volumes. It delivers enhanced clinical flexibility by freeing up time for advanced imaging.

In addition, MAGiC goes beyond the routine, providing complimentary parametric data for a more complete picture.


Go Faster, Further, Deeper with the power of MAGiC.

Faster

With MAGiC you can acquire multiple image contrasts including T1, T2, PD and Inversion Recovery based sequences in just a few minutes. MAGiC also allows you to change image contrasts after acquisition, with simple adjustments – so you never have to worry about acquiring the wrong contrast.




6 scans = 6 contrasts




1 scan = 6 contrasts

Further

With MAGiC dramatically cutting down the time required for a basic scan, you now have the freedom to add on advanced imaging and stay within the same amount of scan time per patient*.



Basic diagnostic acquisitions take over 15 minutes.




Additional acquisitions (ASL, Morphology, DTI and SWAN) acquired in the same amount of time as conventional.

*Results may vary.


Deeper

No longer do you need multiple acquisitions to obtain the parametric data you need. With MAGiC, we ensure B1 uniformity to allow unique results through parametric T1, T2 and PD maps. All of this, at no additional scan time

MAGiC: it's not just an image. It's the full picture.



To acquire parametric MR data, multiple independent sequences are



Within the same single MAGiC scan, complimentary T1, T2, PD, R1 and R2 parametric images are generated.

Faster.
Further.
Deeper.

It's MAGiC!